

## MALE LUER: Integral Lock Ring

### MTLL210 (-1 -3 -4 -5 -6005 -9 -9002 -J2A)

Male Luer Integral Lock Ring to 200 Series Barb, 1/16" (1.6 mm) ID Tubing




### MTLL410 (-9)

Male Luer Integral Lock Ring to 400 Series Barb, 1/16" (1.6 mm) ID Tubing



### MTLL004 (-1 -2 -3 -4 -6005)

Male Luer Integral Lock Ring to 500 Series Barb, 1/16" (1.6 mm) ID Tubing



### MTLL220 (-1 -2 -4 -5 -6005 -9 -9002 -J2A)

Male Luer Integral Lock Ring to 200 Series Barb, 3/32" (2.4 mm) ID Tubing




### MTLL420 (-1 -6005)

Male Luer Integral Lock Ring to 400 Series Barb, 3/32" (2.4 mm) ID Tubing



### MTLL007 (-1 -6005 -9)

Male Luer Integral Lock Ring to 500 Series Barb, 3/32" (2.4 mm) ID Tubing



### MTLL230 (-1 -2 -3 -6005 -9 -9002 -J2A -CM550)

Male Luer Integral Lock Ring to 200 Series Barb, 1/8" (3.2 mm) ID Tubing



### MTLL430 (-3 -6005 -9002 -J2A)

Male Luer Integral Lock Ring to 400 Series Barb, 1/8" (3.2 mm) ID Tubing



### MTLL013 (-1 -2 -3 -6005 -9 -9002)

Male Luer Integral Lock Ring to 500 Series Barb, 1/8" (3.2 mm) ID Tubing



### MTLL240 (-1 -6005 -06 -07 -9 -J2A)

Male Luer Integral Lock Ring to 200 Series Barb, 5/32" (4.0 mm) ID Tubing



### MTLL025 (-1 -2 -4 -5 -6005 -06 -07 -9)

Male Luer Integral Lock Ring to 500 Series Barb, 5/32" (4.0 mm) ID Tubing



### MTLL445 (-6005 -9002)

Male Luer Integral Lock Ring to 400 Series Barb, .170" (4.3 mm) ID Tubing



### MTLL250 (-1 -4 -06 -6005 -07 -8003 -9)

Male Luer Integral Lock Ring to 200 Series Barb, 3/16" (4.8 mm) ID Tubing



### MTLL450 (-9 -9002 -J2A)

Male Luer Integral Lock Ring to 400 Series Barb, 3/16" (4.8 mm) ID Tubing



### MTLL035 (-1 -06 -6005 -07 -9)

Male Luer Integral Lock Ring to 500 Series Barb, 3/16" (4.8 mm) ID Tubing



### MTLL360 (-1 -6005 -9 -9002 -J2A)

Male Luer Integral Lock Ring to 300 Series Barb, 1/4" (6.4 mm) ID Tubing



### MTLL055 (-1 -3 -4 -06 -6005 -9 -9002)

Male Luer Integral Lock Ring to 500 Series Barb, 1/4" (6.4 mm) ID Tubing



## Part Drawings Online!

Nordson MEDICAL provides dimensional drawings available for download for all of its fittings.

[www.nordsonmedical.com](http://www.nordsonmedical.com)

**MALE LUER:**

Slip, or for use with Lock Ring

LUER FITTINGS

**FSLLR (-1 -2 -3 -4 -5 -6005 -07 -9 -9002)**

Finger Snap Luer Lock Ring



**MLRL004 (-1 -6005)**

Male Luer to 500 Series Barb, 1/16" (1.6 mm) ID Tubing (May be used with separate rotating lock ring; FSLLR)



**MLRL007 (-1 -6005 -6012)**

Male Luer to 500 Series Barb, 3/32" (2.4 mm) ID Tubing (May be used with separate rotating lock ring; FSLLR)



**MLSL007 (-9)**

Male Luer to 500 Series Barb, 3/32" (2.4 mm) ID Tubing (May be used with separate stationary lock ring; FSLLR)



**MLRL013 (-1 -6005 -9)**

Male Luer to 500 Series Barb, 1/8" (3.2 mm) ID Tubing (May be used with separate rotating lock ring; FSLLR)



**MLSL013 (-6005)**

Male Luer to 500 Series Barb, 1/8" (3.2 mm) ID Tubing (May be used with separate stationary lock ring; FSLLR)



**MLRL025 (-1 -6005 -9)**

Male Luer to 500 Series Barb, 5/32" (4.0 mm) ID Tubing (May be used with separate rotating lock ring; FSLLR)



**MLRL035 (-1 -6005 -9)**

Male Luer to 500 Series Barb, 3/16" (4.8 mm) ID Tubing (May be used with separate rotating lock ring; FSLLR)



**MLRL055 (-1 -9)**

Male Luer to 500 Series Barb, 1/4" (6.4 mm) ID Tubing (May be used with separate rotating lock ring; FSLLR)



**Need Samples for Evaluation?**

You can call to request most of our standard product in quantities of 5, or less, for free application evaluation.

**+1 (970) 267-5200 Toll Free: (888) 404-5837**

**It's a Lock with Nordson MEDICAL!**

Our MLRL and MLSL Series of fittings can be used both as a slip luer or as a locking luer with the addition of our FSLLR luer lock ring, which allows color coding without flow path contact. Our MLSL/FSLLR combination results in a stationary lock ring. Use a combination of MLRL and FSLLR to create a color coded rotating lock ring.

