


H₂O EXPRESS TEST KIT

FAST, EASY, COST-EFFECTIVE SYSTEM RECOMMENDATIONS

The H₂O Express Test Kit is designed to help you choose the best Pentair® Everpure® solution for your equipment based on the condition of your water!

H₂O Express Test Kit: 180-80193



Testing is easy! Simply complete the instruction form included in the kit. Collect a water sample according to the instructions and send to the address on the form. The sample will be tested and a system recommendation will be sent to the email address(es) provided on the form. Everpure H₂O Express Test Kits are available through an authorized Pentair Dealer.

TEST KIT FEATURES:

- ◆ Evaluation of water quality to determine your treatment needs
- ◆ Mail-in lab test for economical analysis and recommendation

TEST KIT CONTAINS:

- ◆ Test Kit Information Form
- ◆ Three (3) 150 mL Water Sample Vials
 - Please be sure to fill each of the vials full of water, otherwise certain physical characteristics may not be evaluated.
 - To be filled and mailed to Pentair for laboratory testing.
- ◆ Pre-Addressed Padded Return Mailer

TEST RESULTS INCLUDE:

- | | | |
|--------------|-------------|--------------|
| ◆ pH | ◆ TDS | ◆ Hardness |
| ◆ Chlorine | ◆ Silica | ◆ Alkalinity |
| ◆ Chloramine | ◆ Chlorides | ◆ Sulfate |

Test results normally complete in five (5) business days after receipt by Pentair's test lab.

Everpure® H₂O Express Test Customer Information & Test Results Form

The H₂O Test Express Kit contains:

- **Customer Information and Test Results Form.**
- **Three 150mL Water Sample Vials** to be filled and mailed to Pentair for testing. Must include sealable bag to place vials in for shipping.
- **Pre-addressed Mailer** to mail the Customer Information and Test Results Form and Water Sample Vials to Pentair.

INSTRUCTIONS:

- Step 1-** Fill out the section: Test Location, Customer and Dealer Information located below.
- Step 2-** Fill out the section: Food Service Equipment located below describing for what equipment the water is being tested.
- Step 3-** Fill the (3) enclosed vials with cold unfiltered water. Take cold water samples from a spigot at a nearby sink. Run water for three (3) minutes before filling vials.
- Step 4-** Place three (3) bottles in enclosed sealable bag. Enclose the Customer Information and Test Results Form and sealed bag with three (3) sample vials in the pre-addressed envelope.
- Step 5 - Strongly Recommended: Ship by a Trackable Method, i.e. Priority Mail, UPS, or FedEx**

Mail to address labeled on
return envelope:

Pentair H₂O Express Test Lab
101 South Gary Avenue, Suite A
Roselle, IL 60172

Phone: 972.881.9797

WATER TESTING RESULTS

Test Location Name	Address	City	State	Zip Code
Customer Contact Name	Customer Email Address		Customer Phone Number	
Dealer Contact Name	Dealer Email Address		Dealer Phone Number	
Dealer Name	Location	Customer Reference		

TEST LOCATION, CUSTOMER AND DEALER INFORMATION

	Test Results			Test Results	
Total Dissolved Solids (Tds)		ppm		Silica	ppm
Hardness		gpg		pH	units
Chloride		ppm		Total - Chlorine	ppm
Alkalinity		ppm		Free - Chlorine	ppm
Sulfate		ppm		Chloramines (Calculated)	ppm
Pentair Tester Name					Date

FOOD SERVICE EQUIPMENT

RECOMMENDATION NOTES

Equipment Type*	Manufacturer	Model	System	Model	Catalog Page

*List 'Equipment Type' as Combi Oven / Steamer / Ice Machine / Espresso / Coffee Brewer / Etc.

The purpose of this recommendation is to help Pentair users select the best water treatment for their application based upon the water quality parameters presented here and proven technologies to help protect equipment against undesirable water characteristics. Many variables can influence the overall impact water will have on equipment condition or beverage quality including application specifics, equipment design, features, and operator use. Pentair is not liable for equipment condition at any time. Consult with your Pentair Dealer or contact the factory for further assistance selecting and sizing the right Pentair solution for your application.

